

**EEL 4914 Senior Design
Final Design Report**

Fall 2006

Automatic Pill Dispenser

Team: Not a Clue

Submitted by:

B. Bowers – bowersb@ufl.edu – 321.987.7277
M. Preston – mitchelp@ufl.edu – 352.372.5980

I. Project Assumptions and Objectives

Background Information

As the cost of in-home medical care rises, it has become more and more incumbent upon individuals in need of supervised medical care to find a means in which to lower their medical care costs. As such, many individuals who require the administering of many dosages of medications at specific times have turned to devices such as automatic pill dispensers to alleviate the need for an in home nurse on a daily basis. These dispensers range in cost from \$200 up to \$800. These dispensers allow for an in home medical care provider to have a means to regulate a patient's medications without having to constantly supervise the patient. Typical features on these dispensers include automatic pill dispensing at regular intervals, audible warnings, as well as a connection to either a phone line or the internet for monitoring purposes by the medical care provider.

As previously mentioned, many of these devices are rather expensive and can be somewhat cumbersome. Through the use of a simple microprocessor and motor unit an automatic pill dispenser can be produced for a much cheaper price and be much more user friendly.

Project Definition

The team will design and build an automatic pill dispenser. The product will consist of a circular base with 22 fan-like blades that rotate about the central axis. The blades will form the compartments where pills can be manually placed for dispensing at predetermined times. The dispenser will be controlled by a microprocessor that interfaces with an LED display, as well as an alphanumeric keypad that will be utilized as a source for the inputting of data, and selecting from preprogrammed menu items. The user will be able to input the time(s) of day that pills will be dispensed, as well as any warnings and/or precautions that must be followed when the patient takes his or her pills. Lastly, the dispenser will visually and audibly notify the patient when pills are being dispensed, and will also store the time of day that the patient took his or her medications. Finally, the dispenser will automatically adjust the time of the next medication dispersal if necessary, to avoid dosages of medication being taken too closely together. The project will be realized with the development of:

- i. A compact, pill-dispensing unit that can be placed on a table or countertop.
- ii. A microprocessor unit that will control the form and function of the device according to specifications.

Scope of Project

The pill dispenser will be developed with off-the-shelf technology for the design and implementation of the project. The end goal is not to develop any new technologies associated with current manufactured dispensers. Rather, the goal is to design a unit with the same basic functionality, but for a much cheaper price.

Major Objectives

- i. Construct a device that is relatively small and lightweight.
- ii. Develop the software in such a way that patients receive their medication reliably and safely as prescribed by their physician.
- iii. Use as much off the shelf technology, as well as harvest parts from other systems to keep costs low.
- iv. Develop a device that can perform all the necessary functions as stated in the project abstract.

Product Expectations

- i. An audio and visual alarm to notify patient that medication has been dispensed.
- ii. A microprocessor controlled system that will automatically dispense medications at the preset time of day.
- iii. Software that will monitor and record the time that medications have been taken by the patient.
- iv. Software that will automatically adjust future medication dispensation based on when the medication was taken by the patient.
- v. A mechanical locking mechanism that will keep the patient from over medicating.
- vi. An LED display that will give pertinent instructions about the medication to the patient.

II. Customer Requirements

Performance

- Lightweight
- Easy to use for both caretaker and patient
- Well constructed to avoid potential tampering
- Bright warning LED
- 70+ decibel audio warning

Serviceability

- Long life span
- Easy to repair

Features

- Lightweight
- Small dimensions
- Ability to reset the system
- Locking mechanism
- Processor controlled automation for pill dispensing
- Audio and visual alarm
- LED display for patient notification/instructions
- Ability to store times that patient received medications

Reliability

- Stable software
- Mechanical devices encased for safety and durability

Cost

- Low cost
- Off the shelf components that are easily replaced

Safety

- Electrical components encased
- Mechanical parts and motor encased

III. Analysis of Competitive Products

After researching competitive products via the internet, the team made a decision as to what features are necessary to the construction of a useful product. Furthermore, certain features were added that were not typical to current market products. We also assessed the viability of some more complicated features given the amount of time we have to develop out product.

Typical Features of Market Products

- Fully Automatic Pill Dispenser
- Easy Set-up
- Simple to use
- Unlimited # dispenses per day (up to 28 times per day)
- Medication trays
- Lockable with key
- Long Alarm time duration (up to 60 minutes)

Additional Features

There was only one feature in particular that the team felt may not be feasible within the scope of the time given to accomplish the project. Some of the higher end dispensers contain a feature that will notify the patient's caretaker by either phone, e-mail, or both that the patient did not receive their medication within a timely manner. Though this feature is not impossible, it would be quite difficult to reproduce within the time given for design and production of the product.

IV. Concept Selection

In this section, current market product features were taken into account to decide on the features we would include in our pill dispenser. Having taken aforementioned features and concerns into account we decided on the following general concept:

The Programmable Medication Dispenser (PMD) design allows the caregiver to reliably administer medications to a patient without needing to be present every time the medication is scheduled. The PMD allows the caregiver to preprogram up to 21 medication doses through an ergonomically designed interface, utilizing an alphanumeric keypad, and LED display. The basic concept for our project is listed below with a final concept shown in Figure 1.

Hardware Concept Design

Power Supply: The power supply design will provide the necessary power requirements of the PMD. The design requirements are 5 VDC for the microprocessor, and motor controller; additionally, 12 and/or 15 VDC maybe necessary for the motor controller. The current requirements will be mainly dictated by the motor controller design; while the current requirements for the microprocessor are in the mA range, the motor controller may require several amps. Finally, the power supply may require battery backup to avoid loss of user input selections and time keeping functions; depending on the microprocessor and memory designs selected.

Keypad: The keypad input will be a standard 16 key alphanumeric keypad.

Motor: The motor for the PMD will be a stepper motor. The stepper motor selection will be determined by the torque requirements of the dispenser and the number of step divisions required for reliable and accurate medication delivery.

Pill Container and Dispenser: The pill container design will have 22 slots with 21 available for medications. The 21 slots will be labeled so that the caregiver can ensure they are setting the alarm for the correct slot. Additionally, the pill container will incorporate interlock sensors to protect the caregiver and user from harm. Finally, the sensors will provide input to the microprocessor for determining access doors status, medication slot positioning and time feedback of the patient accessing the medications.

Software Concept Design

Microprocessor: The microprocessor will be selected to meet required functionality of the PMD, without wasting money on unneeded features. Additionally, the microprocessor will be chosen so that external memory can be added if necessary, depending on the final algorithm design.


Alarm Module: The alarm module will provide an audible alarm tone. Additionally, the alarm module will provide a visual cue, through the use of LEDs. The design of the alarm module will be to provide a limited alarm function in case of loss of main AC power.

Motor Controller: The design of the motor controller will be determined by the current requirements of the stepper motor selected. The motor controller will take logic inputs from the microprocessor and supply enough current to the stepper motor to meet maximum torque requirements.

LED: The LED display will be a simple black and white, multi-line line display. The LED display will provide information to the caregiver for the purposes of programming and will give feedback on medication compliance by the patient. Additionally, the LED display will convey pre-selected precautions to the patient concerning the medications being currently dispensed.

Complete Design Concept

Programmable Medication Dispenser Block Diagram



V. Project Plans and Scheduling

Our team developed a time line, as well as delegated responsibilities in such a way that all team members must participate in an equal manner. Though we decided to work together and consult as a team on all aspects of the design, certain members were responsible for specific aspects of the project as detailed below.

5.1 **Parts** / Mitchel Preston

All parts were ordered by no later than 9/11/06, and the last part was acquired by 9/20/06.

5.2 **Circuit Design** / Benjamin Bowers

This part of the project was simply comprised of designing a theoretical circuit utilizing software such as PSPICE to test potential circuit implementations for our project. Before anything was soldered, we needed to know that we have a circuit design that will meet the specifications. Furthermore, this aspect of the project helped determine any and all possible challenges met throughout the course of the project given the parts that we decided to utilize.

5.3 **Circuit Implementation** / Benjamin Bowers

This part of the project was the physical implementation of our project once we were satisfied with the results. During this phase of the project, breadboard circuit testing was conducted in order to determine whether our project meets the specifications as well as what was predicted by the software programs. Furthermore, the physical soldering and packaging of the circuit will also take place in this phase of the design once breadboard testing was complete.

5.4 **Product Construction** / Benjamin Bowers (primary) & Mitchel Preston

This phase of the project consisted of the physical building of the unit, itself. Though one member is listed as being primarily responsible, both team members actively participated in the build.

5.5 **Software Implementation** / Mitchel Preston

This phase of the project consisted of the design of the software that controlled the function of the motor, alarm circuit, and memory/LED readout.

Project Timeline (Gantt Chart)

	Wk 1 8/21/06	Wk 2 8/28/06	Wk 3 9/4/06	Wk 4 9/11/06	SpBk 9/18/06	Wk 6 9/19/06	Wk 7 9/20/06	Wk8 9/21/06	Wk9 9/22/06	Wk10 9/23/06	Wk11 9/24/06
Parts /Mitchel P.	X	X	X	X							
Circuit Design / Ben B.	X	X	X	X							
Circuit Implementation /Ben B.				X	X	X					
Product Construction / Ben B. (primary) Mitchel P.						X	X	X	X	X	
Software Implementation / Mitchel P.							X	X	X	X	

VI. Unit Cost

One of the goals of this design was to produce a product comparable to current market pill dispensers but at a more affordable cost. The lowest cost unit that was found through internet research was approximately \$150. Though cheaper units can be bought, the aforementioned unit was the cheapest one that the team could locate with features similar to our design. The cost for the team's design can be seen in the chart below. However, it must be noted that many of the components were provided by Mike Stapleton, free of charge.

Design Unit Cost Analysis Table

Item	Cost per unit
Current Limiting Diodes	Free from Mike
Stepping Motor	\$10.00
Resistors	\$3.00
L293 Motor Driver	\$1.00
Plastiboard	\$5.00
3/4 inch Dowel Rod	\$0.50
555 Timer	\$1.69
PNP transistors	\$2.59
Bright Red LED	\$2.69
75 dB Buzzer	\$3.29
Capacitors	Free from Mike
4 MHz Crystal	Free from Mike
D1307 Chip	Free from Mike
32 Crystal	Free from Mike
Lithium Battery	Free from Mike
Battery Holder	Free from Mike
PIC16F877A I/P	\$7.50
Switching Diode N4914	Free from Mike
Cookie Tin	\$2.00
Screws	\$6.00
LCD	Free from Mike
16 Key Alphanumeric Keypad	MicroP
Total Unit Cost:	\$45.26

The total unit cost was well under \$150. However, this number does not take into account the cost of many minor items such as batteries, capacitors, and crystal chips. Even so, the total unit cost would still be well under the original goal. As such, the team's design is an extremely affordable unit.

VII. Completed Design Analysis

7.1 Design Changes

There were very few changes to the original design. Originally, our team wanted to be able to program the module for all 21 slots. However, it was later discovered that the PIC16F877 I/P chip contained only 8kB of memory. As such, there was not enough memory to allow for this feature after all of the other software design features were implemented. Therefore, the team decided to allow for the ability to program two slot times in order to show proof of concept understanding and implementation. The chart below shows lists all of the original design features with a verification mark next to those that were met.


Feature	Functional	Not Functional
Ability to Keep Time	X	
Audiovisual Alarm	X	
Ability to Set Multiple Alarm Times	X	
View Past Alarms	X	
View Current Alarm Times	X	
Auto Positioning of the Carousel	X	
Time Battery Backup	X	
Backup Alarm Memory	X	
Backup Alarm Adjust	X	
Keypad Programmability	X	
LCD Functional	X	

7.2 Design Problems/Bugs


Despite taking more time than was expected, we ran into very few major problems or bugs during the process of designing and building this unit. For the physical construction of the unit, the only problem that arose involved the carousel. The carousel rests on a vertical, wooden rod that is connected to the motor gear. As such, when the processor steps the motor gear the rod is rotated 15° which, in turn, also rotates the carousel to the next pill slot. However, the team discovered a problem with the rod not being perfectly vertical. This caused some swiveling of the carousel that resulted in the auto-positioning sensor to not always be activated properly. However, this problem was easily resolved with some manual bending of the can to compensate for the imperfection on the rod's positioning. On the software side, the limited amount of memory within the PIC chip forced the team to truncate some of the original design parameters. The team had originally hoped to produce a unit that could program all 21 pill slots. However, the iteration in the code that was required to achieve this was very taxing on the memory space. Therefore, the team decided that instead of cutting some of the other design features, it would be prudent to simply allow for the programming of a few slots to demonstrate the concept. This decision allowed for the continued development of all other features to maintain robustness within the unit design.


7.3 Software Logic Flowcharts

The program for the medication dispenser consisted of a loop with three main routines.




The Timer Related Block Logic is as seen below:






The Sensor Block Logic is as follows:



Input User Block Logic:



VIII. Appendix

8.1 PMD Code

Cntr VAR BYTE EEProm	'This holds the index that will be used to read out of
-------------------------	--

Cntr=1

CntrIndex VAR BYTE	'This holds the index to control the number of reads at beginning of program stored at \$01
--------------------	---

'AlmPntr stored at \$00

```
EEPromIndex VAR BYTE    'This holds the EEPROM storage address
```

EEPromIndex = \$02

VAR BYTE

Minutes VAR BYTE[21] 'The actual time is memory spot 0, alarm 1st slot is memory spot 1 with the remaining 20 following.

Hours VAR BYTE[21]

Date VAR BYTE[21]

Month VAR BYTE[21]

Year VAR BYTE[21]

```
Prcn VAR BYTE[21]
```

'Array for tracking precaution to use as messsage

AlmFlag VAR BIT[21]

'Flag for tracking if alarm has been set and if alarm

has gone off

'This array may not work because I can't

initialize the values and if power is lost

'any initialization routine would undo the set

flags. I can initialize the values when the alarm slots

'are set and that may work out?

AlmOn VAR BIT

AlmSetMin VAR BYTE[21]

'Holds the minutes that the alarm was set for

AlmSetHour VAR BYTE[21]

' Holds the hours that the alarm was set for

MinMedsTaken VAR BYTE [21]

'Memory to record minutes when patient took meds

HrMedsTaken VAR BYTE[21]

'Memory to record hours when patient took meds

TempHours VAR BYTE 'This variable is used by the hours setting routine to strip the 12/24 and AM/PM bits off.

'Additionally, it will be used to hold the

difference between the current time meds taken and

'the next alarm time.

TimeHours	VAR BYTE	This variable is used to retain the hours once the 12/24 and the AM/PM bits are stripped off
-----------	----------	--

CurrentAlm VAR BYTE 'This will hold the current alarm for comparison
 NextAlarm VAR BYTE 'This will hold the next alarm for comparison
 CurrentAlmAmPm VAR BIT 'For indicating if the time is AM or PM for adjusting the time between alarms.
 NextAlmAmPm VAR BIT 'For indicating if the time is AM or PM for adjusting the time between alarms.
 CurrentAlmTens VAR BIT 'For indicating if the time has a 10 hour bit this is needed because of the registers format
 NextAlmTens VAR BIT 'For indicating if the time has a 10 hour bit this is needed because of the registers format
 DiffBtwnHrs VAR BYTE 'This will hold the difference between the two alarm times for adjusting the second alarm time
 CurrentAlmDay VAR BYTE 'This will hold the current alarm's day for comparison
 NextAlmDay VAR BYTE 'This will hold the next alarm's day for comparison
 NewAlmTime VAR BYTE 'This will hold the new alarm time
 DiffDay VAR BIT

ReviewSlot VAR BYTE 'Variable for reviewing when patient took meds by pointing to the memory spot being viewed
 MinTimeBtwnMeds VAR BYTE 'This variable holds the selected amount of times between meds

TempNum VAR BYTE 'This variable holds and returns the 2 digit number from the Keys subroutine
 TempNum2 VAR BYTE[2] 'For the temp storage of entered 2 digit numbers
 SecDig VAR BIT
 Digit VAR BYTE 'This variable is used by the Keys subroutine to track
 'the number of digits entered TO

control the Return function
 Key VAR BYTE
 TempAlmIndx VAR BYTE 'This variable is used by next slot to track slot being set to update the TimeAlmIndx which is reset
 'by the menu each time it is entered
 TimeAlmIndx VAR BYTE 'This displays the slot number of the current slot being programmed see the day, hour, etc.
 AlmPntr VAR BYTE 'This points at the current alarm waiting to go off

SettingAlm VAR BIT 'This variable indicates if the alarm is being set
 DisplayPrctn VAR BYTE 'This bit is used to determine which precaution to display when subroutine is called

PORTB.5=1
PORTC.2=1 'Alarm module line

'Once the code is fully written this code can be deleted

'I2CWrite PORTC.4, PORTC.3, \$d0, \$00, [Seconds] 'Setting the CH bit7 to 0
ensures oscillator operation.
'Pause 50

'I2CWrite PORTC.4, PORTC.3, \$d0, \$01, [Minutes[0]]
'Pause 50

'I2CWrite PORTC.4, PORTC.3, \$d0, \$02, [Hours[0]]
'Pause 50

'I2CWrite PORTC.4, PORTC.3, \$d0, \$04, [Date[0]]
'Pause 50

'I2CWrite PORTC.4, PORTC.3, \$d0, \$05, [Month[0]]
'Pause 50

'I2CWrite PORTC.4, PORTC.3, \$d0, \$06, [Year[0]]
'Pause 50

'End of deletion, the code can be deleted because we don't want the value in the clock chip being reinitialized with each power up.

'Without this section when the unit is repowered the clock will need to be set the first time but everytime after that the time should
'be retained by the D1307 chip.

PORTD.0=1 'This section initializes the values for the keypad
PORTD.1=1
PORTD.2=1
PORTD.3=1

AlmOn=0

Read \$00,AlmPntr
Read \$01,CntrIndex
StartLoop: IF (CntrIndex>=Cntr) AND (EEPromIndex<\$FF) Then
Read EEPROMIndex, Prectn[Cntr]
EEPromIndex=EEPromIndex+1
Read EEPROMIndex,Minutes[Cntr]
EEPromIndex=EEPromIndex+1
Read EEPROMIndex, Hours[Cntr]

```
EEPromIndex=EEPromIndex+1
Read EEPROMIndex,Date[Cntr]
EEPromIndex=EEPromIndex+1
Read EEPROMIndex,Month[Cntr]
EEPromIndex=EEPromIndex+1
Read EEPROMIndex,Year[Cntr]
EEPromIndex=EEPromIndex+1
Cntr=Cntr+1
GoTo StartLoop
ENDIF
```

'Increments the counter'

```
MainLoop:  I2CRead PORTC.4, PORTC.3, $d0, $00, [Seconds]
          I2CRead PORTC.4, PORTC.3, $d0, $01, [Minutes[0]]
          I2CRead PORTC.4, PORTC.3, $d0, $02, [Hours[0]]
          I2CRead PORTC.4, PORTC.3, $d0, $04, [Date[0]]
          I2CRead PORTC.4, PORTC.3, $d0, $05, [Month[0]]
          I2CRead PORTC.4, PORTC.3, $d0, $06, [Year[0]]
          TimeHours=Hours[0]&%00011111
```

```

'IF alarm1=0 Then           'This executes as long as the alarm is not
going off
    IF (Hours[0]>=$40)AND(Hours[0]<=$52) Then LCDOut $FE,1,"Time:
AM ",HEX2 TimeHours,":",HEX2 Minutes[0],":", HEX2 Seconds

```

```
IF(Hours[0]>=$60)AND (Hours[0]<=$72) Then LCDOut $FE,1,"Time:
PM ",HEX2 TimeHours,":",HEX2 Minutes[0],":", HEX2 Seconds
```

```
LCDOut $FE,$C0,HEX2 Month[0],"/",HEX2 Date[0],"/",HEX2 Year[0]
'EndIF
```

Pause 300 'To cutdown on LCD flash

```
***** Alarm Test and Set Off
*****
*****
```

IF (Minutes[0]==Minutes[AlmPntr])
AND(Hours[0]==Hours[AlmPntr])AND(Date[0]==Date[AlmPntr])AND(Month[0]==Month[AlmPntr])AND (Year[0]==Year[AlmPntr])AND (alarm1=0)AND(AlmOn=0) Then

IF (PORTC.0=1) AND (PORTC.1=0) Then

```

DisplayPrctn=Prctn[AlmPntr]      'This sets the precaution
display variable equal to the one chosen for this alarm
GoSub Precautions                  'This jumps to the
precautions subroutine and displays the chosen precaution
AlmSetMin [AlmPntr]=Minutes[AlmPntr] 'This records the
minutes of the current alarm time
AlmSetHour[AlmPntr]=Hours[AlmPntr] 'This record
the minutes of the current alarm time

AlmOn=1

GoSub MtrDvr

PORTC.2=0      'This port turns on the alarm module
alarm1=1       'Temp variable to turn off alarm for test will
use Alarm Flag later
EndIF
EndIF

***** Alarm turn off routine
*****
IF (PORTC.1=1)AND(AlmOn=1) Then      'This executes if the
pill output door was opened turning off the alarm
PORTC.2=1      'Turns off alarm

I2CRead PORTC.4, PORTC.3, $d0, $01, [MinMedsTaken[AlmPntr]]
'Records the minutes when the patient took the meds
I2CRead PORTC.4, PORTC.3, $d0, $02, [HrMedsTaken[AlmPntr]]
'Records the hours the patient took the meds

Minutes[AlmPntr]=Minutes[AlmPntr]-1 'This subtracts one minute from
the current alarm to prevent alarm from going off again
'for
this same alarm time

AlmOn=0      'This resets to prevent this loop from being
entered if the alarm is not going off

LCDOut $FE,1,"Alarm off"
Pause 2000
GoSub Precautions
Pause 2000

*****Alarm adjustment
section*****

```

```

'           CurrentAlmAmPm=HrMedsTaken[AlmPntr]&%00100000 'This is now. If
this equals 1 then the time was PM and I need to add 12 hours to the TimeHours
'           CurrentAlmTens=HrMedsTaken[AlmPntr]&%00010000 'This is now. If
this equals 1 then the hours has a tens component
'           CurrentAlm=HrMedsTaken[AlmPntr]&%00001111 'This gets the current
alarm time actual one's hours
'           CurrentAlmDay=Date[AlmPntr]

```

```

'***** This code puts the current
alarm time in 24 hour format *****

```

```

'12PM in 24 hour format
'           IF
(CurrentAlmAmPm=%00100000)AND(CurrentAlmTens=%00010000)AND(CurrentAl
m=2)Then CurrentAlm=CurrentAlm+%00001010 'add 10
'10PM to 11PM in 24 hour format
'           IF
(CurrentAlmAmPm=%00100000)AND(CurrentAlmTens=%00010000)Then
CurrentAlm=CurrentAlm+%000010110 'add 22
'1PM to 9PM in 24 hour format
'           IF (CurrentAlmAmPm=%00100000)Then
CurrentAlm=CurrentAlm+%00001100 'add 12
'12AM in 24 hour format
'           IF
(CurrentAlmAmPm=%00000000)AND(CurrentAlmTens=%00010000)AND(CurrentAl
m=2)Then CurrentAlm=0
'10AM to 11AM in 24 hour format
'           IF
(CurrentAlmAmPm=%00000000)AND(CurrentAlmTens=%00010000)AND(CurrentAl
m=0)Then
'           CurrentAlm=10
'           Else
'           CurrentAlm=11
'           EndIF
'1AM to 9AM in 24 hour format require no actions

```

```

'           NextAlmAmPm=Hours[AlmPntr+1]&%00100000 'This is next. If this
equals 1 then the time was PM and I need to add 12 hours to the TimeHours
'           NextAlmTens=Hours[AlmPntr+1]&%00010000 'This is next. If this
equals 1 then the hours has a tens component
'           NextAlarm=Hours[AlmPntr+1]&%00001111 'This gets the next
alarm time actual one's hours

```

```

'           NextAlmDay=Date[AlmPntr+1]

'           ***** This code puts the next
alarm time in 24 hour format *****
'           '12PM in 24 hour format
'           IF
(NextAlmAmPm=%00100000)AND(NextAlmTens=%00010000)AND(NextAlarm=2)T
hen NextAlarm=NextAlarm+%00001010 'add 10
'           '10PM to 11PM in 24 hour format
'           IF (NextAlmAmPm=%00100000)AND(NextAlmTens=%00010000)Then
NextAlarm=NextAlarm+%00010110 'add 22
'           '1PM to 9PM in 24 hour format
'           IF (NextAlmAmPm=%00100000)Then
NextAlarm=NextAlarm+%00001100 'add 12
'           '12AM in 24 hour format
'           IF
(NextAlmAmPm=%00000000)AND(NextAlmTens=%00010000)AND(NextAlarm=2)T
hen NextAlarm=0
'           '10AM to 11AM in 24 hour format
'           IF
(NextAlmAmPm=%00000000)AND(NextAlmTens=%00010000)AND(NextAlarm=0)T
hen
'           NextAlarm=10
'           Else
'           NextAlarm=11
'           EndIF
'           '1AM to 9AM in 24 hour format require no actions

'           DiffDay=CurrentAlmDay-NextAlmDay           'results = 0 if the days
are the same

'           IF (DiffDay<>0)AND(NextAlmAmPm=%00000000)Then
NextAlarm=NextAlarm+24 'The next alarm is the following day and AM => +24
'           'This makes the next equation workout

'           DiffBtwnHrs=NextAlarm-CurrentAlm           'This gets the difference
between the current alarm and the next alarm so
'           'I can now figure out if the next alarm needs adjustment

'           IF (MinTimeBtwnMeds>=DiffBtwnHrs) Then
DiffBtwnHrs=MinTimeBtwnMeds-DiffBtwnHrs 'DiffBtwnHrs now has the adjustment
margin

'           IF DiffBtwnHrs<>0 Then LCDOut $FE,1, "Delay Meds by ", HEX2
DiffBtwnHrs, " Hrs"

```

***** Including alarm adjustment code to even this point cause
 problems with power up display and causes the alarms
 ***** malfunction in memory this code takes the code to \$1E88
 which is very close to the end of the available
 ***** memory therefore I suspect this code is overwriting other code
 such as the stack
 ***** without extra memory we will not be able to include the alarm
 adjustment feature.

' IF DiffBtwnHrs=0 Then GoTo AlmOK 'No adjustment to the next
 alarm is necessary

***** If the days are different (DiffDay not =0) and the times are both
 AM or PM. The time span is >12 hrs and ***

***** no alarm
 adjustment is necessary therefore no actions required ****

' IF
 (DiffDay<>0)AND(CurrentAlmAmPm=%00100000)AND(NextAlmAmPm=%00100000
) Then GoTo AlmOK

' IF
 (DiffDay<>0)AND(CurrentAlmAmPm=%00000000)AND(NextAlmAmPm=%00000000
) Then GoTo AlmOK

' NewAlmTime=CurrentAlm+DiffBtwnHrs
 'Future = Now +DiffBtwnHrs

***** Current and Next alarm Same day Current
 alarm AM New alarm time 0-11 (AM) *****

' IF (DiffDay=0)&(CurrentAlm<=11)&(NewAlmTime<=11) Then GoSub
 AlmSub

'This step is not necessary because can't have a zero unless zero was added
 which would have caused the code to skip this

'IF NewAlmTime=0 Then NewAlmTime=12 'This takes the 24
 hour clock and converts to 12 hour clock

'IF NewAlmTime>9 Then
 'NewAlmTime=NewAlmTime-10 'This strips of the 10 therefore the
 result is 0 or 1
 'NewAlmTime=NewAlmTime|%01010000 'This is now 12hr,AM,10hr,0
 or 1

'EndIF

'NewAlmTime=NewAlmTime|%01000000 'This is now 12hr,AM,1-9

```

'EndIF

***** Current and Next alarm Same day Current
alarm AM/PM New alarm time 12-23 (PM) *****
'      IF (DiffDay=0)AND(NewAlmTime>=12)AND(NewAlmTime<=23) Then
'      '&(CurrentAlm<=11)
'      '      NewAlmTime=NewAlmTime-12
'      '      IF NewAlmTime=0 Then NewAlmTime=12 'This takes the 24 hour clock
and converts to 12 hour clock

'      '      IF NewAlmTime>9 Then
'      '      NewAlmTime=NewAlmTime-10 'This strips of the 10 therefore the
result is 0,1 or 2
'      '      NewAlmTime=NewAlmTime|%01110000 'This is now 12hr,PM,10hr,0 or
1

'      EndIF

'      NewAlmTime=NewAlmTime|%01100000 'This is now 12hr,PM,1-9
'      EndIF

***** Current and Next alarm Same day Current
alarm PM New alarm time 12-23 (PM) *****
'      IF (DiffDay=0)&(NewAlmTime>=12)&(NewAlmTime<=23) Then
'      '&(CurrentAlm>11)
'      '      NewAlmTime=NewAlmTime-12 'This puts the 24 hour clock back to
12 hour format
'      '      IF NewAlmTime=0 Then NewAlmTime=12 'This takes the 24
hour clock and converts to 12 hour clock

'      '      IF NewAlmTime>9 Then
'      '      NewAlmTime=NewAlmTime-10 'This strips of the 10 therefore the
result is 0,1 or 2
'      '      NewAlmTime=NewAlmTime|%01110000 'This is now 12hr,PM,10hr,0
or 1

'      EndIF

'      NewAlmTime=NewAlmTime|%01100000 'This is now 12hr,PM,1-9

'      EndIF

***** Current and Next alarm Same day Current
alarm PM New alarm time 0-11 (AM next day) *****
'      IF (DiffDay=0)AND(NewAlmTime>=24)AND(NewAlmTime<=35) Then
'      '&(CurrentAlm>11)
'      '      NewAlmTime=NewAlmTime-24

```



```

'          Date[AlmPntr+1]=Date[AlmPntr+1]+1          'This should
increment the date one day
'          GoSub AlmSub
          'IF NewAlmTime=0 Then NewAlmTime=12          'This takes the 24
hour clock and converts to 12 hour clock

          'IF NewAlmTime>9 Then
          'NewAlmTime=NewAlmTime-10    'This strips of the 10 therefore the
result is 0,1 or 2
          'NewAlmTime=NewAlmTime|%01010000 'This is now 12hr,AM,10hr,0
or 1

          'EndIF

          'NewAlmTime=NewAlmTime|%01000000 'This is now 12hr,AM,1-9
'          EndIF

          ***** Current and Next alarm Diff days Current
alarm PM New alarm time 0-11 (AM next day) *****
          'Compiler indicated that the processor ran out of memory with this section
of code above

'AlmSub: IF NewAlmTime=0 Then NewAlmTime=12    'This takes the 24 hour clock
and converts to 12 hour clock

'          IF NewAlmTime>9 Then
'          NewAlmTime=NewAlmTime-10    'This strips of the 10 therefore the
result is 0,1 or 2
'          NewAlmTime=NewAlmTime|%01010000 'This is now 12hr,AM,10hr,0
or 1

'          EndIF

'          NewAlmTime=NewAlmTime|%01000000 'This is now 12hr,AM,1-9
'          Return

          *****
*****
'AlmOK:    pause 10

          AlmPntr=AlmPntr+1
          IF AlmPntr=22 Then AlmPntr=1          'Prevents the alarm pointer
from going out of bounds

```

' IF AlmPntr<TimeAlmIndx Then AlmPntr=AlmPntr+1 'This
increments the alarm pointer to point at the next alarm in line

' IF DiffBtwnHrs<>0 Then

' Hours[AlmPntr] = NewAlmTime
' LCDOut \$FE,1,"Hours next med",HEX2 Hours[AlmPntr]
' Pause 3000
' EndIF

alarm1=0

EndIF

***** End of Alarm Module

' IF MinMedsTaken[AlmPntr-1]<>Minutes[0] Then alarm1=0 'This will
temporarily allow the alarm to function for the second alarm

'it will have to be replaced by the alarm flag.This routine has a
problem

'with the last memory slot because the AlmPntr is not advanced
therefore

'it is comparing with the memory slot from previous alarm which
will always

'make the if condition true

'IF AlmOn=1 Then
'GoSub Precautions
'Pause 2000
'EndIF

'Loop1:

IF PORTC.0=0 Then
LCDOut \$FE,1,"Close Main Lid"
Pause 2000
'Switch Debounce time delay
'IF PORTC.0=0 Then GoTo Loop1
EndIF

'Loop2:

IF (PORTC.1=1) AND (AlmOn=0) Then
LCDOut \$FE,1,"Close Output Door"

```

        Pause 2000
        'Switch Debounce time delay
    ' IF PORTC.1=1 Then GoTo Loop2
    EndIF

'Loop3:   IF PORTC.6=0 Then LCDOut $FE,1,"Close Pill "
'         Pause 10
'         'Switch Debounce time delay
'         IF PORTC.6=0 Then GoTo Loop3

'Loop4:   IF PORTC.7=0 Then LCDOut $FE,1,"At first slot"
'         'Pause 10
'         'Switch Debounce time delay
'         IF PORTC.7=0 Then GoTo Loop4

MenuIndex=0

PORTD.3=0
IF PORTD.4=0 Then
    Pause 100                                'Switch debounce
    PORTD.3=1
    GoTo Menu                                'This enters the menu loop
EndIF
PORTD.3=1                                'This turns off the column because menu was not
selected

PORTD.2=0
IF PORTD.7=0 Then 'Pound key pushed because last alarm was set
    Pause 100                                'Switch debounce
    PORTD.2=1
PillsRdy: GoSub MtrDvr                        'This will execute until slot1 is positioned to be the
next slot over hole
    IF PORTC.7=0 Then
        LCDOut $FE,1,"All alarms set"
        Pause 2000
        GoTo MainLoop
    Else
        GoTo PillsRdy
    EndIF
EndIF
PORTD.2=1                                'This turns off the column because the subroutine
was not selected

Write $00,AlmPntr

```

GoTo MainLoop

'1: Pause 10

Menu: LCDOut \$FE,1,"Set time? Enter"

 'LCDOut \$FE,\$C0,"Next, Back, or Escape"

 Pause 1000 '1 sec delay to prevent last loop keypress from
making selections in this loop

 MenuIndex=1 'Initializes the menu choice subroutine's index
number

 'GoSub ChooseEnter

 PORTD.0=0

 IF PORTD.7=0 Then 'Enter was pressed

 Pause 200

 PORTD.0=1

 TimeAlmIndx=0

 PORTD.0=1

 GoTo Minute

 EndIF

 PORTD.0=1

 GoSub MenuChoice

 'PORTD.3=0

 'IF PORTD.4=0 Then GoTo Menu

 'IF PORTD.5=0 Then GoTo MedsTaken 'Back was pressed

 'IF PORTD.6=0 Then GoTo AlSlot1 'Next was pressed

 'IF PORTD.7=0 Then GoTo MainLoop 'Escape was pressed

 'PORTD.3=1

 GoTo Menu

'2: Pause 10

AlSlot1:LCDOut \$FE,1,"Set 1st slot alarm. Enter"

 'LCDOut \$FE,\$C0,"Next, Back, or Escape"

 DisplayPrctn=0

 Pause 1000

 'GoSub ChooseEnter 'Enter was pressed

 PORTD.0=0

 IF PORTD.7=0 Then

 Pause 200

```

PORTD.0=1
CntrIndex=1
Write $01,CntrIndex
EEPromIndex=$02 'Resets the EEPROMIndex
AlmPntr=1 'This points to the first alarm slot in memory
TimeAlmIndx=1 'This points to the 1st alarm memory slots
TempAlmIndx=1 'This variable is being set so that it will be
incremented correctly each time next slot is entered.

```

SettingAlm=1 'This variable tells the minute, hour, day, month, and year
not to write the values entered to the D1307 chip

```

MtrLoop:GoSub MtrDvr
IF PORTC.6<>0 Then GoTo MtrLoop
LCDOut $FE,1,"At first slot "
Pause 2000
IF SettingAlm=1 Then GoSub PrectnChoice

```

GoTo Minute

EndIF 'Ends the if statement for choosing to set the

1st alarm slot

```
PORTD.0=1
```

```

GoSub MenuChoice
'PORTD.3=0 'These would execute if the customer doesn't want to set the
first slot alarm

```

```
'IF PORTD.4=0 Then GoTo Menu
```

```
'IF PORTD.5=0 Then GoTo Menu 'Back was pressed
```

```
'IF PORTD.6=0 Then GoTo NextAlm 'Next was pressed
```

```
'IF PORTD.7=0 Then GoTo MainLoop 'Escape was pressed
```

```
'PORTD.3=1
```

GoTo AlSlot1

```

'3:      Pause 10
NextAlm:LCDOut $FE,1,"Set next alarm"
Pause 1000

```

```
PORTD.0=0
```

```
IF PORTD.7=0 Then 'Enter was chosen
```

```
Pause 200
```

```

PORTD.0=1          'Turns off first column if enter was chosen
SettingAlm=1 'This variable tells the minute, hour, day, month, and year
not to write the values entered to the D1307 chip
TempAlmIndx=TempAlmIndx+1 'This increments the temp alarm
index each time the next slot is selected
LCDOut $FE,1,"Setting Alm Slot",HEX2 TempAlmIndx
'Pause 2000
'IF TempAlmIndx=22 Then          'Out of slots therefore this ends the
option of setting further times
'LCDOut $FE,1,"All slots have been set"
'Pause 2000
'GoTo MainLoop
'ENDIF          'End out of slots if statement

TimeAlmIndx=TempAlmIndx
CntrIndex=TimeAlmIndx          'This advances the counter index for
the read cycle on power up
Write $01,CntrIndex
GoSub MtrDvr          'This advances the carousel one slot
GoSub PrecnChoice

GoTo Minute
ENDIF          'Ends the enter if then statement

PORTD.0=1          'Turns off the first column if enter was not chosen
GoSub MenuChoice
'PORTD.3=0 'These would execute if the customer doesn't want to set the
first slot alarm
'IF PORTD.4=0 Then GoTo Menu

'IF PORTD.5=0 Then GoTo AlSlot1 'Back was pressed
'IF PORTD.6=0 Then GoTo TimeBtwnMeds          'Next was pressed
'IF PORTD.7=0 Then GoTo MainLoop          'Escape was pressed
'PORTD.3=1
GoTo NextAlm

'4:          Pause 10
TimeBtwnMeds:LCDOut $FE,1,"Set minimum time"
LCDOut $Fe,$C0,"between meds"
Pause 1000
PORTD.0=0 'Turns on the first column of keys
IF PORTD.7=0 Then 'Enter was pressed

PORTD.0=1          'Turns off the enter key
Pause 100

```

```

LCDOut $FE,1,"Enter min number of"
LCDOut $Fe,$C0,"hours 0-12 btwn doses"
GoSub Keys
MinTimeBtwnMeds=TempNum
LCDOut $FE,1,"Hours chosen is",HEX2 MinTimeBtwnMeds
Pause 3000 'Gives time to read the
message
IF (MinTimeBtwnMeds<0)|(MinTimeBtwnMeds>7)Then GoTo
TimeBtwnMeds 'Entry was out of bounds
GoTo MainLoop
EndIF
PORTD.0=1 'Turns off the enter key because enter wasn't
chosen
GoSub MenuChoice
'PORTD.3=0 'These would execute if the customer doesn't want
to set the first slot alarm
'IF PORTD.4=0 Then GoTo Menu

'IF PORTD.5=0 Then GoTo NextAlm 'Back was pressed
'IF PORTD.6=0 Then GoTo ReviewAlmTime 'Next was
pressed
'IF PORTD.7=0 Then GoTo MainLoop 'Escape was pressed
'PORTD.3=1
GoTo TimeBtwnMeds

'5: Pause 10
ReviewAlmTime: LCDOut $FE,1,"View Alarm Times"
LCDOut $Fe,$C0,"Enter,Next,Back,Esc"
Pause 1000
PORTD.0=0 'Turns on the first column of keys
IF PORTD.7=0 Then 'Enter was pressed

PORTD.0=1
Pause 1000
LCDOut $FE,1,"Enter slot to review"

GoSub Keys
ReviewSlot=TempNum
IF ReviewSlot>21 Then GoTo ReviewAlmTime

TimeHours=Hours[ReviewSlot]&%00011111

```

```
IF (Hours[ReviewSlot]>=$40)AND(Hours[ReviewSlot]<=$52)
Then LCDOut $FE,1,HEX2 ReviewSlot,"Slot Set: AM ",HEX2 TimeHours,":",HEX2
Minutes[ReviewSlot]
```

```
IF(Hours[ReviewSlot]>=$60)AND (Hours[ReviewSlot]<=$72)
Then LCDOut $FE,1,HEX2 ReviewSlot,"Slot Set: PM ",HEX2 TimeHours,":",HEX2
Minutes[ReviewSlot]
```

```
LCDOut $FE,$C0,HEX2 Month[ReviewSlot],"/",HEX2
Date[ReviewSlot],"/",HEX2 Year[ReviewSlot]
```

```
Pause 3000
```

```
GoTo MainLoop
```

```
EndIF
```

```
PORTD.0=1 'Turns off the enter key if enter was not pressed
```

```
GoSub MenuChoice
```

```
'PORTD.3=0 'Turns on the 4th column of keys
```

```
'IF PORTD.5=0 Then GoTo TimeBtwnMeds 'Back was
pressed
```

```
'IF PORTD.6=0 Then GoTo MedsTaken 'Next
was pressed
```

```
'IF PORTD.7=0 Then GoTo MainLoop
```

```
'Escape was pressed
```

```
'PORTD.3=1
```

```
GoTo ReviewAlmTime
```

```
'6: Pause 10
```

```
MedsTaken:LCDOut $FE,1,"View time meds taken"
```

```
Pause 1000
```

```
PORTD.0=0 'Turns on the first column of keys
```

```
IF PORTD.7=0 Then 'Enter was pressed
```

```
PORTD.0=1
```

```
Pause 1000
```

```
LCDOut $FE,1,"Enter slot to review"
```

```
GoSub Keys
```

```
ReviewSlot=TempNum
```

```
'IF ReviewSlot>21 Then GoTo MedsTaken
```

```
TimeHours=AlmSetHour[ReviewSlot]&%00011111
```



```
IF
(AlmSetHour[ReviewSlot]>=$40)AND(AlmSetHour[ReviewSlot]<=$52) Then LCDOut
$FE,1,HEX2 ReviewSlot,"Set: AM ",HEX2 TimeHours,":",HEX2
AlmSetMin[ReviewSlot]
```

```
IF(AlmSetHour[ReviewSlot]>=$60)AND
(AlmSetHour[ReviewSlot]<=$72) Then LCDOut $FE,1,HEX2 ReviewSlot,"Set: PM
",HEX2 TimeHours,":",HEX2 AlmSetMin[ReviewSlot]
```

```
TimeHours=HrMedsTaken[ReviewSlot]&%00011111
```

```
IF
(HrMedsTaken[ReviewSlot]>=$40)AND(HrMedsTaken[ReviewSlot]<=$52) Then
LCDOut $Fe,$C0,HEX2 ReviewSlot,"Taken: AM ",HEX2 TimeHours,":",HEX2
MinMedsTaken[ReviewSlot]
```

```
IF(HrMedsTaken[ReviewSlot]>=$60)AND
(HrMedsTaken[ReviewSlot]<=$72) Then LCDOut $Fe,$C0,HEX2 ReviewSlot,"Taken:
PM ",HEX2 TimeHours,":",HEX2 MinMedsTaken[ReviewSlot]
```

```
Pause 3000
```

```
GoTo MainLoop
```

```
EndIF
```

```
PORTD.0=1 'Turns off the enter key if enter was not pressed
```

```
' GoTo ReviewLoop 'This loop executes until user presses
escape to return to Main Loop
```

```
GoSub MenuChoice
```

```
'PORTD.3=0 'These would execute if the customer doesn't want to set the
first slot alarm
```

```
'IF PORTD.4=0 Then GoTo Menu
```

```
'IF PORTD.5=0 Then GoTo ReviewAlmTime
```

```
'Back was pressed
```

```
'IF PORTD.6=0 Then GoTo Menu
```

```
'Next was
```

```
pressed
```

```
'IF PORTD.7=0 Then GoTo MainLoop
```

```
'Escape was
```

```
pressed
```

```
'PORTD.3=1
```

```
GoTo MedsTaken
```

```
*****
*****
```

'Below this point are all the subroutines: Time and Alarm setting, Keys, Precautions, Menu Choices

```
Minute:    LCDOut $FE,1,"Enter minute's"
           'Pause 2000
           'LCDOut $FE,1,HEX2 TimeAlmIndx
           Pause 1000
           GoSub Keys                'Calls Keys subroutine after subroutine
program returns here
           Minutes[TimeAlmIndx]=TempNum
           IF SettingAlm=0 Then I2CWrite PORTC.4, PORTC.3, $d0, $01,
[Minutes[TimeAlmIndx]]
           'The SettingAlm variable tells this part of the code if the alarm is being set
           'if the alarm is being set then this is skipped to avoid setting the clock with
the alarm time.
```

```
SubLoop1:LCDOut $FE,1,"Is",HEX2 Minutes[TimeAlmIndx],"Correct?"
           LCDOut $Fe,$C0,"Enter=Yes and 4=No"
           Pause 200
```

```
           PORTD.0=0
           IF PORTD.7=0 Then 'Enter was pressed for yes
           Pause 200
           PORTD.0=1
```

```
           'LCDOut $FE,1,Minutes[0]
           'LCDOut $Fe,$C0,Minutes[1]
```

```
           'Pause 1000
           IF SettingAlm=1 Then
           Write EEPromIndex, Minutes[TimeAlmIndx]
           EEPromIndex=EEPromIndex+1
           EndIF
           GoTo Hour
EndIF
```

```
           IF PORTD.5=0 Then
           Pause 200
           GoTo Minute
EndIF
```

```
           GoTo SubLoop1
```

```
Hour:    LCDOut $FE,1,"Enter hour's"
           Pause 200
```

```

        GoSub Keys                                'Calls Keys subroutine after subroutine
program returns here
        Hours[TimeAlmIndx]=TempNum

SubLoop2:LCDOut $FE,1,"Is",HEX2 Hours[TimeAlmIndx],"Correct?"
        LCDOut $Fe,$C0,"Enter=Yes and 4=No"
        Pause 200

        PORTD.0=0
        IF PORTD.7=0 Then 'Enter was pressed for yes
            Pause 200

            GoTo Subloop3

        EndIF

        IF PORTD.5=0 Then
            Pause 200
            GoTo Hour
        EndIF

        GoTo SubLoop2

Subloop3: LCDOut $FE,1,"AM=1 and PM=4"
        Pause 200
        IF PORTD.4=0 Then 'AM was pressed
            Pause 200
            TempHours=$40
            LCDOut $FE,1,"AM was selected"
            Pause 1000
            GoTo Subloop4
        EndIF

        IF PORTD.5=0 Then
            Pause 200
            TempHours=$60
            LCDOut $FE,1,"PM was selected"
            Pause 1000
            GoTo Subloop4
        EndIF
        GoTo Subloop3

Subloop4: PORTD.0=1                                'Disables the first column of numbers
        Hours[TimeAlmIndx]=TempHours|Hours[TimeAlmIndx]
        IF SettingAlm=0 Then I2CWrite PORTC.4, PORTC.3, $d0, $02,
[Hours[TimeAlmIndx]]

```

'The SettingAlm variable tells this part of the code if the alarm is being set
'if the alarm is being set then this is skipped to avoid setting the clock with
the alarm time.

```
IF SettingAlm=1 Then
    Write EEPromIndex, Hours[TimeAlmIndx]
    EEPromIndex=EEPromIndex+1
EndIF
```

Day: LCDOut \$FE,1,"Enter Day"

Pause 200

GoSub Keys

'Calls Keys subroutine after

subroutine program returns here

Date[TimeAlmIndx]=TempNum

IF SettingAlm=0 Then I2CWrite PORTC.4, PORTC.3, \$d0, \$04,
[Date[TimeAlmIndx]]

'The SettingAlm variable tells this part of the code if the alarm is being set
'if the alarm is being set then this is skipped to avoid setting the clock with
the alarm time.

SubLoop5:LCDOut \$FE,1,"Is",HEX2 Date[TimeAlmIndx],"Correct?"

LCDOut \$Fe,\$C0,"Enter=Yes and 4=No"

Pause 200

'GoSub ChooseEnter

PORTD.0=0

IF PORTD.7=0 Then 'Enter was pressed for yes

Pause 200

PORTD.0=1

IF SettingAlm=1 Then

Write EEPromIndex, Date[TimeAlmIndx]

EEPromIndex=EEPromIndex+1

EndIF

GoTo Mnth

EndIF

IF PORTD.5=0 Then

Pause 200

GoTo Day

EndIF

GoTo SubLoop5

Mnth: LCDOut \$FE,1,"Enter Month"

```

        Pause 200
        GoSub Keys                                'Calls Keys subroutine after
subroutine program returns here
        Month[TimeAlmIndx]=TempNum
        IF SettingAlm=0 Then I2CWrite PORTC.4, PORTC.3, $d0, $05,
[Month[TimeAlmIndx]]
        'The SettingAlm variable tells this part of the code if the alarm is being set
        'if the alarm is being set then this is skipped to avoid setting the clock with
the alarm time.

```

```

SubLoop6:LCDOut $FE,1,"Is",HEX2 Month[TimeAlmIndx],"Correct?"
        LCDOut $Fe,$C0,"Enter=Yes and 4=No"
        Pause 200

```

```

        PORTD.0=0
        IF PORTD.7=0 Then 'Enter was pressed for yes
                Pause 200
                PORTD.0=1

```

```

                IF SettingAlm=1 Then
                Write EEPromIndex, Month[TimeAlmIndx]
                EEPromIndex=EEPromIndex+1
                EndIF

```

```

                GoTo Yr
        EndIF

```

```

        IF PORTD.5=0 Then
                Pause 200
                GoTo Mnth
        EndIF
        GoTo SubLoop6

```

```

Yr: LCDOut $FE,1,"Enter Year"
        Pause 200
        GoSub Keys                                'Calls Keys subroutine after subroutine
program returns here
        Year[TimeAlmIndx]=TempNum
        IF SettingAlm=0 Then I2CWrite PORTC.4, PORTC.3, $d0, $06,
[Year[TimeAlmIndx]]
        'The SettingAlm variable tells this part of the code if the alarm is being set
        'if the alarm is being set then this is skipped to avoid setting the clock with
the alarm time.

```

```

        SubLoop7:LCDOut $FE,1,"Is",HEX2 Year[TimeAlmIndx],"Correct?"

```

```
LCDOut $Fe,$C0,"Enter=Yes and 4=No"  
Pause 200
```

```
PORTD.0=0  
IF PORTD.7=0 Then 'Enter was pressed for yes  
Pause 200  
PORTD.0=1
```

```
IF SettingAlm=1 Then  
Write EEPromIndex, Year[TimeAlmIndx]  
EEPromIndex=EEPromIndex+1  
EndIF
```

```
SettingAlm=0 'This resets the alarm set to false so the time  
or alarm can be changed in the future.  
GoTo MainLoop  
EndIF
```

```
IF PORTD.5=0 Then  
Pause 200  
GoTo Yr  
EndIF  
GoTo SubLoop7
```

```
*****  
*****
```

Keys: PORTD.0=0

```
IF PORTD.4=0 Then  
Pause 200 'switch depressed debounce  
TempNum2[Digit]=$01  
Key=$01  
EndIF
```

```
IF PORTD.5=0 Then  
Pause 200 'switch depressed debounce  
TempNum2[Digit]=$04  
Key=$04 'switch release debounce
```

```
EndIF
```

```
IF PORTD.6=0 Then  
Pause 200 'switch depressed debounce
```

```

TempNum2[Digit]=$07
Key=$07                                'switch release debounce
EndIF

'IF PORTD.7=0 Then Key="*"
IF (PORTD.4=0)OR(PORTD.5=0)OR(PORTD.6=0)Then
    KLoop1:        LCDOut $FE,1,HEX Key
    Pause 200      'give the lcd time to
update
    IF (PORTD.4=0)OR(PORTD.5=0)OR(PORTD.6=0) Then GoTo
KLoop1
    Digit=Digit+1
    Pause 200      'switch release
debounce
EndIF
PORTD.0=1

PORTD.1=0
IF PORTD.4=0 Then
    Pause 200      'switch depressed debounce
    TempNum2[Digit]=$02
    Key=$02
EndIF
IF PORTD.5=0 Then
    Pause 200      'switch depressed debounce
    TempNum2[Digit]=$05
    Key=$05
EndIF

IF PORTD.6=0 Then
    Pause 200      'switch depressed debounce
    TempNum2[Digit]=$08
    Key=$08
EndIF
IF PORTD.7=0 Then
    Pause 200      'switch depressed debounce
    TempNum2[Digit]=$00
    Key=$00
EndIF
IF
(PORTD.4=0)OR(PORTD.5=0)OR(PORTD.6=0)OR(PORTD.7=0) Then
    KLoop2:        LCDOut $FE,1,HEX Key

```

```

                                Pause 200                                'give the lcd time to
update
                                IF
(PORTD.4=0)OR(PORTD.5=0)OR(PORTD.6=0)OR(PORTD.7=0) Then GoTo KLoop2

                                Digit=Digit+1
                                Pause 200      'switch release debounce
                                EndIF

PORTD.1=1

PORTD.2=0
IF PORTD.4=0 Then
    Pause 200                                'switch depressed debounce
    TempNum2[Digit]=$03
    Key=$03
EndIF
IF PORTD.5=0 Then
    Pause 200                                'switch depressed debounce
    TempNum2[Digit]=$06
    Key=$06
EndIF
IF PORTD.6=0 Then
    Pause 200                                'switch depressed debounce
    TempNum2[Digit]=$09
    Key=$09
EndIF
'IF PORTD.7=0 Then LCDOut $FE,1,"#"
IF (PORTD.4=0)OR(PORTD.5=0)OR(PORTD.6=0) Then
KLoop3:      LCDOut $FE,1,HEX Key
    Pause 200                                'give the lcd time to
update
                                IF (PORTD.4=0)OR(PORTD.5=0)OR(PORTD.6=0) Then GoTo
KLoop3
                                Digit=Digit+1
                                Pause 200
                                EndIF                                'switch release debounce
PORTD.2=1

'    PORTD.3=0
'    IF PORTD.4=0 Then LCDOut $FE,1,"Menu"
'    IF PORTD.5=0 Then LCDOut $FE,1,"Back"
'    IF PORTD.6=0 Then LCDOut $FE,1,"Next"
'    IF PORTD.7=0 Then LCDOut $FE,1,"Escape"

```



```

'    PORTD.3=1

'    Pause 10

    IF (Digit=0)OR(Digit=1) Then
        GoTo Keys
    Else
        Digit =0
    EndIF

    TempNum= TempNum2[0]
    TempNum=TempNum<<4           'Shifts the low byte of TempNum to
the High byte of TempNum
    TempNum=TempNum|TempNum2[1]   'ORs the High byte of
temp num with the low byte of 1

    Return

```

MenuChoice: PORTD.3=0 'These would execute if the customer doesn't want to set the first slot alarm

```

    IF PORTD.5=0 Then           'Back was pressed
        MenuIndex=MenuIndex-1
        IF MenuIndex=<0 Then MenuIndex=6
    EndIF

    IF PORTD.7=0 Then   'Escape was pressed
        PORTD.3=1
        GoTo MainLoop
    EndIF

    IF PORTD.6=0 Then           'Next was pressed
        MenuIndex=MenuIndex+1
        IF MenuIndex=>7 Then MenuIndex=1
    EndIF

    PORTD.3=1

    IF MenuIndex=1 Then GoTo Menu
    IF MenuIndex=2 Then GoTo AlSlot1
    IF MenuIndex=3 Then GoTo NextAlm
    IF MenuIndex=4 Then GoTo TimeBtwnMeds
    IF MenuIndex=5 Then GoTo ReviewAlmTime
    IF MenuIndex=6 Then GoTo MedsTaken

```

Return

Precautions: Pause 100 'delay to keep lcd from flashing

IF DisplayPrctn=0 Then LCDOut \$FE,1,"No Precautions"

IF DisplayPrctn=1 Then
LCDOut \$FE,1,"Causes Drowsiness"
LCDOut \$Fe,\$C0,"Use no alcohol"
EndIF

IF DisplayPrctn=2 Then LCDOut \$FE,1,"Take with water"

IF DisplayPrctn=3 Then LCDOut \$FE,1,"Take with food"

IF DisplayPrctn=4 Then
LCDOut \$FE,1,"Don't take with"
LCDOut \$Fe,\$C0,"nitrates."
EndIF

IF DisplayPrctn=5 Then LCDOut \$FE,1,"May cause dizziness"

IF DisplayPrctn=6 Then
LCDOut \$FE,1,"Don't use with"
LCDOut \$Fe,\$C0,"herbal products"
EndIF

IF DisplayPrctn=7 Then
LCDOut \$FE,1,"check before using"
LCDOut \$Fe,\$C0,"with OTC Meds"
EndIF

Return

PrectnChoice:Pause 100

Subloop8: LCDOut \$FE,1,"Choose precaution"

Pause 2000 'Time to read message
GoSub Precautions

```

Pause 2000 'Gives time to read precaution
LCDOut $FE,1,"Is that correct"

Pause 2000 'Gives time to read message
PORTD.0=0 'Turns on first column of keys

pressed IF PORTD.7=0 Then 'Precaution was chosen and enter was
pressed
Pause 100 'Debounce pause
Prectn[TimeAlmIdx]=DisplayPrctn 'Enters the chosen precaution
message in the precaution array
PORTD.0=1 'Turns off first column of keys

IF SettingAlm=1 Then
Write EEPromIndex, Prectn[TimeAlmIdx]
EEPromIndex=EEPromIndex+1
EndIF

Return 'The precaution has been chosen and this returns the
program to the gosub call in alarm set
'GoTo Minute
EndIF

selection PORTD.3=0 'Turns on last column of keys for precaution

'IF PORTD.4=0 Then GoTo Menu

IF PORTD.5=0 Then DisplayPrctn=DisplayPrctn+1 'Next was pressed
'Loops the DisplayPrctn counter around to the previous precaution

'EndIF
IF PORTD.6=0 Then DisplayPrctn=DisplayPrctn-1 'Back was pressed
'Loops the DisplayPrctn counter around to the next precaution

'EndIF
'IF PORTD.7=0 Then GoTo MainLoop 'Escape was pressed
PORTD.3=1 'Turns off the last column of keys
IF (DisplayPrctn<0)Then DisplayPrctn=7
IF (DisplayPrctn>7) Then DisplayPrctn=0
GoTo Subloop8

```

MtrDvr: MtrIndex=MtrIndex+1 'This increments the index causing it to toggle each time this subroutine is called

```
IF MtrIndex=1 Then
PORTB.0=1'Chip enable clockwise
PORTB.2=1
PORTB.1=0 'Motor advance
Pause 1000 'Pause to allow carousel to settle
PORTB.1=1 'Turn off motor
PORTB.0=0 'Turn off enable
EndIF
```

```
IF MtrIndex=0 Then
PORTB.0=1'Chip enable clockwise
PORTB.2=0
PORTB.1=1 'Motor advance
Pause 1000 'Pause to allow carousel to settle
PORTB.2=1 'Turn off motor
PORTB.0=0 'Turn off enable
EndIF
```

Return